



ASPIRE *Cristalle*

The fight against breast cancer is
not one-dimensional

FUJIFILM
Value from Innovation



3D Mammography System featuring an ergonomic design for comfortable patient positioning with a 4-second tomo sweep. The Cristalle comes with several features to improve the technologist experience such as, one-touch automated mirrored positioning, one-touch exam start and One-Shot Phantom M Plus that allows for weekly QC to be completed in as little as seven minutes.

Patient Comfort The patented Comfort Paddle gently and evenly contours to the breast shape to provide a noticeably more comfortable exam with more precise compression. Comfort Comp is a software feature that automatically readjusts compression to minimize discomfort. Stationary face shield, patient calming mood lighting and accent decals to create better patient experiences are standard.

iAEC combines image analysis technology with intelligent sampling to automatically optimize techniques and processing for individual tissue characteristics and breast densities, streamlining workflow while improving image quality.

Detector features an innovative Hexagonal Close Pattern (HCP) pixel design with low noise, high speed switching technology engineered to deliver higher acquisition efficiency for finer detail, higher DQE, and higher MTF. Easy start up and sleep modes extend detector life up to 300%.