



Persona C

Full size C-arm - go anywhere design with 20 kW of power

FUJIFILM
Value from Innovation

Mobile fluoroscopy solution with removable grid and dedicated 'radiography mode' for high quality x-ray imaging.

Designed for easy, precise navigation and imaging across a broad range of diagnostic orthopedic imaging and minimally invasive surgical procedures and pain management.

Automated low dose organ programs. DSA digital subtraction, vascular package for angiography

Product Highlights

- Streamlined design with light easy versatility for tight workspaces
- Optimized for surgical guidance
- Surgeon's console with Dual 21" displays
- dRF detector 30 or 21cm, large FOV, low dose and high image quality
- 3 acquisition modes: low dose, high quality & snapshot
- 3 mode Active Cooling lowers tube heat as much as 50% (option)

Category exclusive features

- Dual lasers at tube & detector
- Stay in place arm positioning
- Easy lift front wheels
- 81cm largest opening in its class
- Only system w/+65° to -95° orbital rotation
- Lightest, up to 200 lbs. lighter
- Only battery free system

