



Persona C

*Full Size C-Arm - go anywhere,
do anything with 20 kW of power*

Mobile fluoroscopy solution with removable grid and dedicated 'radiography mode' for high quality still imaging. Designed for easy, precise navigation and imaging across a broad range of diagnostic imaging and minimally invasive surgical procedures.

- Full-size high-power multi-function
- Optimized for surgical guidance
- Separate tech and surgeon's consoles
- Dedicated dual display surgical viewing console
- Built on DR flat panel technology for ultra-low dose fluoroscopy and vascular imaging
- Streamlined design with light easy versatility for tight work spaces



Product Highlights

- dRF detector 30 cm or 21 cm, large FOV, low dose and high image quality
- Dedicated radiography mode high resolution static x-rays, simplify procedure time pre- & post- imaging with one system
- Light, easy positioning, large 81 cm arc for easy access and precise positioning
- Automated low dose organ programs
- DSA digital subtraction, vascular package for angiography

FUJIFILM Healthcare Americas Corporation

81 Hartwell Avenue, Suite 300, Lexington, MA 02421

fujifilmhealthcare.com

800.431.1850

© FUJIFILM Healthcare Americas Corporation

DOC-0057429-A

FUJIFILM
Value from Innovation