



SCENARIA View

*A powerful premium performance
128-slice CT solution*

Scenaria View CT is for radiology departments, ED's and outpatient centers that need dependable radiology applications capabilities with available advanced applications modules for interventional CT, extended coverage Shuttle Scanning (for perfusion exams), cardiac CTA and dual energy examinations.



Scenaria View features a novel lateral shift table with ± 10 cm shift range and Fujifilm's advanced iterative reconstruction, Intelli IPV for dose reduction with natural image texture.

Product Highlights

- Table capacity: 550 lbs.
- Bore Size: 80 cm
- Slice Count: 64/128
- Field of View: 50 cm
- Rotation Speed: 0.35 sec
- Metal Artifact Reduction: HiMAR Plus (standard)

FUJIFILM Healthcare Americas Corporation

81 Hartwell Avenue, Suite 300, Lexington, MA 02421

fujifilmhealthcare.com

800.431.1850

© FUJIFILM Healthcare Americas Corporation

DOC-0057429-A

FUJIFILM
Value from Innovation