



THE
IMAGING
SOLUTIONS COMPANY™

AMRAD  Medical

PARTNERS IN IMAGING EXCELLENCE

DFMT ELITE

Premier
Radiographic
System

Where
Simplicity
Meets
Functionality



- DFMT Elite is fully featured and can be configured with a variety of generators and X-ray tubes to meet the clinical requirements of any imaging environment
- Traditional floor mount design means no learning curve for the technologist
- Industry leading Warranty, Customer Service and Technical Support
- The DFMT Elite is compatible with all DR solutions

Designed for high volume throughput for hospitals, orthopedics and medical imaging facilities. The DFMT Elite can perform a broad range of exams while optimizing technologist productivity and patient workflow. Simplicity and functionality provide broad positioning capability for patient procedures, including cross table laterals. Feature-rich controls and ease of use enables the technologist to focus on patient care.

Performance, Versatility & Value



Overhead Patient Handgrip

For lateral exams and weight bearing studies for 360° positioning with increments at every 5° for maximum flexibility

Remote Collimator Light

Conveniently turn collimator light on from the wall stand

Smart Button Technology

Rotating the angulation handle from table to wall stand, the buttons for vertical and horizontal movement intuitively swap functions for seamless operation

Table Height Indicator

Ensures optimal working height

Adaptive Optimal Display

As the angulation handle is rotated to the wall stand, the display changes orientation for easy viewing

Simple Service

LED operational status indicators

T.E.C.H.

Table Elevation Control Handle for table top lock release and vertical movements

Summit Industries & MXR Imaging

Through our Partners in Imaging Excellence Program, MXR aligns its portfolio of services with manufacturing leaders in medical imaging. In partnership with Summit Industries and the AMRAD® Medical brand, a recognized leader providing high-quality equipment with an exceptional return on investment, our companies go beyond just equipment by providing seamless coordination, exceptional service and the earliest access to new products, technologies and upgrades. Our consultative approach to working with diagnostic imaging centers and urgent care facilities is unmatched. Innovation and customer satisfaction are at the core of everything we do.



4909 Murphy Canyon Road, Suite 120
San Diego, CA 92123
Phone: 888-278-9933
www.mxr imaging.com

