



Persona CS

Compact all-in-one C-arm for ASCs, orthopedic & pain management

Advanced C-arm solution designed for rapid single unit positioning even within the tightest OR environments. With 2.8-5 kW power this compact c-arm is perfectly balanced to provide exceptional imaging performance without compromising maneuverability.

Automated low dose organ programs. DSA digital subtraction, roadmapping, vascular package for angiography



FUJIFILM
Value from Innovation

Product Highlights

- Streamlined design with light easy versatility for tight workspaces
- Compact all-in-one system with integrated surgical viewing & technologist control displays
- Largest opening and significantly lighter than typical systems in its class
- 30 cm or 21 cm dRF detector, large FOV, low dose, high image quality
- Dedicated radiography mode produces high resolution static x-rays for pre- & post- imaging to simplify procedure times with one system
- Light, easy positioning, large 81 cm opening for easy access & precise positioning

Category exclusive features

- Dual lasers at tube & detector
- Stay in place arm positioning
- 81 cm largest opening in its class
- Only system w/+65° to -95° orbital rotation
- Only battery free system