



Persona CS

*Compact all-in-one c-arm for ASCs,
orthopedic & pain management*

Advanced c-arm solution designed for rapid single unit positioning within even the tightest OR environments. 5 kW of power this compact c-arm is perfectly balanced to provide superior x-ray performance without compromising on maneuverability.

- Compact all-in-one design with integrated surgical viewing & tech control displays
- Larger opening and significantly lighter than typical systems in its class.
- Streamlined design with light easy versatility for tight work spaces



Product Highlights

- dRF detector 30 cm or 21 cm, large FOV, low dose high IQ
- dRF dedicated radiography mode high resolution static x-rays, simplify procedure time pre- & post- imaging with one system
- Light, easy positioning, large 81 cm arc for easy access & precise positioning
- Automated Low Dose Organ Programs
- DSA digital subtraction, vascular package for angiography

FUJIFILM Healthcare Americas Corporation

81 Hartwell Avenue, Suite 300, Lexington, MA 02421

fujifilmhealthcare.com

800.431.1850

© FUJIFILM Healthcare Americas Corporation

DOC-0057429-A

FUJIFILM
Value from Innovation