

MAX Series C-Arm Tables



Easy c-arm access

500 lb weight capacity

2" patient comfort pad

Light weight and easy to move

24" wide top

Made in the USA

Standard Accessories:

Hand Control



Patient Straps




Optional Accessories



Standard Side Rails



Paper Roll



Foot Control

*Streamline 1, 2, & 3 require upgrade to use both hand and foot controls simultaneously



Modular Side Rails



Surgical Arm Boards

*Require side rails



Procedural Arm Boards

Configurations/Specifications

MAX-ERCP Series

MAX4-ERCP

Height 30"-40"
Lateral Tilt 17°
Trendelenburg 13°
Longitudinal Travel 16"


MAX5-ERCP

Height 33.5-43.5"
Lateral Tilt 17°
Trendelenburg 13°
Longitudinal Travel 16"
Lateral Travel 10"

Standard MAX Series


MAX 1

Height 28"-38"




MAX 2

Height 28"-38"




MAX 3

Height 28"-38"



MAX 4

Height 30"-40"



MAX 5

Height 33.5"-43.5"

