

VMAX Series C-Arm Tables



Instant-Response 4-way float

Easy c-arm access

500 lb weight capacity

2" patient comfort pad

Light weight and easy to move

24" wide top

Made in the USA

Standard Accessories:

Joystick



Patient Straps



Fixed Side Rails



Optional Accessories



Procedural Arm Boards



Paper Roll



Foot Control



Catheter Tray



Modular Side Rails



Surgical Arm Boards
*Require siderails



Rectangular Arm Board



Lollipop Arm Board



*Trendelenberg only available for VMAXT-Pro

Configurations/Specifications



VMAX-Pro

Height 28"-38"
4 Way Float
Longitudinal Travel 36" (+/-) 18"
Lateral Travel 14" (+/- 7")

VMAXT-Pro

Height 30"-40"
4 Way Float
Longitudinal Travel 36" (+/-) 18"
Lateral Travel 14" (+/- 7")
Trendelenberg +/- 15°

