



# OASIS **VELOCITY**

## AI Accelerates High performance Open MRI

**FUJIFILM**  
Value from Innovation

Experience the patient satisfaction advantages of Deep Learning Reconstruction (DLR – Rise) combined with true high field Open MRI. Lightning fast imaging and exquisite image quality now combine with unmatched comfort.

Velocity's Deep Convolutional Neuro Network reconstruction engine, workflow engineered RF coils, open sided architecture and unique lateral table movement bring valuable benefits to your technologists, radiologists and patients.

Differentiate your practice, increase your volume, enhance your imaging excellence with Velocity.

### Product Highlights

- **DLR Rise: Deep Learning Reconstruction**  
-intelligent noise reduction
- **IP-RAPID acquisition with enhanced SNR for exceptional image quality**
- **Head/Spine Integrated RF coil technology**
- **16 channel RF system**
- **33 mT/m, 100 T/m/s gradients**
- **Blanket coil technology provides maximum convenience, flexibility, comfort and performance**
- **82 cm wide table with 20 cm lateral movement and 660 lbs capacity accommodates patients of all sizes**

