



VETSMART

The Digital Imaging System That Thinks For You



2020imaging.net/vetsmart



VETSMART


The latest in Veterinary X-Ray technology



- Fixed 40" SID
- Toshiba x-ray tube 140 HU 1-2mm focal spot
- kVp range: 40-125kV
- Single Phase 208-230VAC | Minimum 30 amp service
- High frequency 32 kW capacitor assisted 400mA generator
- 15" Touch Screen Operator Console



- Exclusive table base/table top designs
 - 40" base w/ 64" standard table top
 - 40" base w/ 48" or 56" table top options
 - Oval top w/urine track
- Unique Push and Open pull out tray to access DR detector
- Exclusive Push and Open Front panel latch mechanism for easy generator access



- Low input amperage
- SMART APR
- Remote real-time monitoring, calibration and maintenance features (optional)
- Exclusive self-diagnostic tool based on lot and machine learning
- Eliminate need for onsite service technicians



**ELIMINATE HAVING TO
MEASURE THE PATIENT**

- Patented ATD (Auto Thickness Detection) technology
- Open APR integration: control kV, mA, and mS
- Exclusively control generator dosage (embedded AI/lot) by detecting body part for ideal exposure
- Automatic calculation of patient thickness

Veterinary Specific Software

ACQUISITION INTERFACE



- Under 6 second ⌚ Processing
- Background Supression
- Preset X-ray Exams
- Animal Specific
- Anatomical Program (APR)
- Fully Customizable
- Advanced Image Processing
- Image Adjustment Tools
- Touchscreen Enabled
- User Friendly Interface

Multi-funtion Studylist Advanced Viewer

- PM Software Integration*
- Burn/Import CDs
- Custom Search/Query
- Customized User Settings
- Vertibral Heart Score Tool incl.
- Fast, User Friendly Interface
- Image Fine Tuning
- Exam Room Viewing



MOBILE VIEWING WITH M-OPAL[†]

iPad®, iPhone®, Android™, Blackberry®, any mobile or smart device!



- Lightweight & Fast
- Patient List & Viewer
- Simple Annotations Enabled
- Tap, Drag, and Pinch
- Patient Filter & Search

17x17 UNIVERSAL DR FLAT PANEL



- Tethered Interface
- Cesium Scintillator
- 16-bit Images
- Cassette Sized
- 139 micron | 3.6 lp/mm
- 5 year MFG warranty



Advanced
Superior
Efficient
Reliable
Affordable

2020imaging.net/vetsmart

