



THE
IMAGING
SOLUTIONS COMPANY™



KONICA MINOLTA

PARTNERS IN IMAGING EXCELLENCE

Compact Portable X-Ray System

An Imaging
System that
Moves with You




Improve Your Workflow

The specially modified version of ULTRA image acquisition software improves workflow by presenting vital decision-making information for the operator, helping facilitate faster, more confident acquisition.

Features

A full array of detector holders and grids are available to protect and enhance the performance of your mobile or portable X-Ray.

The Konica Minolta Portable X-Ray system is compact, lightweight and fully equipped to take radiology on the road. The compact design is agile and lightweight (100 lbs.) The 3.0kW generator offers fully regulated mA/kVp and is suitable for core portable radiology applications. The built-in laptop and cassette storage allows the operator to operate in tight spaces.



ULTRA + DR Portable

The ULTRA+ DR portable kit features the lightweight wireless cesium panel along with the reliable SR 130 portable X-Ray system. The digital package uses a mounted laptop tray and convenient panel slide to hold the laptop and the panel.

Service Coverage

MXR Imaging provides National Service Support for multiple modalities. This includes coverage for CT, PET/CT, MRI, Ultrasound, PACS/RIS, and General Radiology Imaging equipment. We offer customizable service offerings to healthcare networks, hospitals, imaging centers and Independent Service Organizations (ISOs), as well as Orthopedic, Neurology, Cardiology and Veterinary clinics.

ULTRA + Wireless DR Specifications

Weight: 7.9 lbs (CSI)
Carbon Fiber and Mg Alloy Construction
Pixel Matrix-2304 x 2800
3.7 lp/mm (CSI)
Cassette-sized
Lithium Ion Battery (2 Included)

SR 130 X-Ray System Specifications

3.0kw Output-Fully Regulated mA/kVp
KVP Range 40-100
MA Range 15 and 30
MAS Range 0.15-120



Konica Minolta & MXR Imaging

Together, MXR Imaging and Konica Minolta go beyond just imaging. These companies stand out by providing seamless coordination, exceptional service and the earliest access to new products, technologies and upgrades. MXR and Konica Minolta focus on innovation and customer satisfaction at the core of everything they do.



4909 Murphy Canyon Road, Suite 120
San Diego, CA 92123
Phone: 888-278-9933
www.mxrimaging.com

