



Barrier Technologies[®]
Your Safety is Our Concern[™]

Radiation Protective Eyewear



2026

Your Safety is our Concern™

Welcome to our 2025 Catalog

As you can see, we have revamped our catalog to make it easier for you to use. We have instituted a new RX reference system which will help you easily identify which RX types can be accommodated by each frame.



Plano



Single Vision




BiFocal



Progressive

Through our dedication to innovation, we have made a significant investment in Free-Form Robotic glass surfacing technology which allows us to provide anti-fatigue Single Vision, high definition Bi-Focal, and personalized Progressive prescription with up to a +4.50 add power.



Our partnership with Corning and their [Med-X®](#) line of 0.75mm Pb protective lenses enables us to offer you, our customer, the clearest vision possible while providing maximum protection. By adding the options for ClearView AR (anti-reflective) and AF (anti-fog) coatings to all of our available styles, our offerings in designs and materials cannot be surpassed.



We invite you explore our catalog to see these innovations and how we continue to set the standard for quality and customer service for Radiation Protective Leaded Eyewear.

Ray-Ban®

∞ Meta

NEW



The Most Advanced Radiation Protection Glasses Available

Ray-Ban Meta glasses combine iconic style with the latest in wearable tech, allowing you to capture photos and videos, listen to music with the built-in Bluetooth speakers, interact with AI, and stay connected with calls or texts using just voice commands.

CAPTURE

Capture your unique perspective with the new ultra-wide 12MP camera and five-mic system.

LISTEN

Discreet open-ear speakers let you enjoy rich audio quality, music and calls, while keeping you connected to the world around you.

CALL

Make calls and send messages on your smartphone, Whatsapp, Messenger and SMS, completely hands-free – simply by using voice commands.
+ POV/Videocall : a unique way to call and share your point of view handsfree.

AI*

Spark your creativity, get information and control your glasses with your voice - powered by Meta AI.

With Multimodal AI, you can use your glasses to look and ask Meta AI to take a photo and give information based on what you are looking at to identify objects, translate and caption.

LIVESTREAM

Share your unique view of special moments by livestreaming them through your glasses to Instagram and Facebook.

Prescription Lenses



Meta Ray-Ban AI glasses are prescription lens compatible.

Ray-Ban® Eyewear Collection

Ray-Ban®



BALORAMA

- Weight: 74 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Colors: Black

Ray-Ban®



JACKIE OHH

- Weight: 69 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Colors: Black

Ray-Ban®



JUSTIN

- Weight: 65 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Colors: Black

Ray-Ban®



NEW WAYFARER

- Weight: 59 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Colors: Black

Oakley® Leaded Eyewear



GASCAN

- O-Matter™ stress-resistant frame material is both lightweight & durable for all-day comfort and protection
- Three-Point Fit holds lenses in precise optical alignment
- Weight: 74 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Polished Black, Matte Black, Multi-Camo



FUELCELL

- Unobtainium™ sure-grip inset zones provide durability & performance
- O-Matter™ stress-resistant frame material is both lightweight & durable for all-day comfort and protection
- Three-Point Fit holds lenses in precise optical alignment
- Weight: 72 Grams
- Frame Colors: Polished Black, Multi-Camo



TWOFACE

- Unique pairing of lightweight O-Matter™ upper & high strength Aluminum lower rims
- Lens positioning offers superior comfort, clarity, & razor-sharp vision at every angle
- O-Matter™ stress-resistant frame material is both lightweight & durable for all-day comfort and protection
- Weight: 64 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Polished Black w/Steel Grey, Triple Black



FIVE SQUARED

- O-Matter™ stress-resistant frame material is both lightweight & durable for all-day comfort and protection
- Three-Point Fit holds lenses in precise optical alignment while eliminating pressure points
- Weight: 58 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Colors: Polished Black

Oakley® Leaded Eyewear



CROSSLINK XL

- Sleek and modern look matches Oakley performance innovation with a comfortable fit
- Constructed of light weight, durable O-Matter™ material
- Unobtainium™ nosepads and ear socks for enhanced grip
- Weight: 62 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Satin Black, Grey



JUPITER SQUARED

- True Oakley attitude with refinement
- Unique Oakley Riveting on the temple arm
- O-Matter™ frame ensures comfort during extended procedures
- Unobtainium® nose piece & earsocks for enhanced grip
- Weight: 65 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Colors: Polished Black



MAINLINK

- Lightweight O-Matter® frame for optimal comfort
- Three-Point fit holds the lenses for precise optical alignment
- Inspired by the Crosslink for enhanced Rxability
- Weight: 66 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Colors: Matte Black



HOLBROOK RX

- O-Matter™ stress-resistant frame material is both lightweight & durable for all-day comfort and protection
- Three-Point Fit holds lenses in precise optical alignment
- Weight: 65 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Satin Black, Polished Clear

COSTA Leaded Eyewear



BRINE



- Durable lightweight TR-90 nylon frame construction
- Hydrolite hypoallergenic rubber prevents slippage
- Vented ports promote maximum airflow and reduce fogging
- Fits small to medium faces
- Weight: 56 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Colors: Black, Tortoise



BLACKFIN



- Durable lightweight TR-90 nylon frame construction
- Hydrolite hypoallergenic rubber prevents slippage
- Fits medium to large faces
- Weight: 73 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Colors: Matte Black



FANTAIL



- Modified (DBL) eliminates frame resting on face
- Hydrolite hypoallergenic rubber prevents slippage
- Fits medium to large faces
- Weight: 69 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Colors: Black

COSTA® Leaded Eyewear




DREAM

- Ideal for all Rx types
- Fits large to extra large faces
- Weight: 69 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Colors: Black



PERMIT

- Modified (DBL) eliminates frame resting on face
- Hydrolite hypoallergenic rubber prevents slippage
- Vented ports promote maximum airflow and reduce fogging
- Fits medium to large faces
- Weight: 69 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Colors: Matte Black



TUNA ALLEY

- Modified (DBL) eliminates frame resting on face
- Hydrolite hypoallergenic rubber prevents slippage
- Vented ports promote maximum airflow and reduce fogging
- Fits medium to large faces
- Weight: 69 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Colors: Matte Black, Tortoise



Leaded Eyewear



P-17

- Peak Performance Ballistic Eyewear
- ANSI Z87.1 rated frames
- Innovative design for Unparalleled Protection
- Designed for Superior Protection & Performance
- Weight: 56 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Gloss Black



GRAVITY



- Peak Performance Ballistic Eyewear
- ANSI Z87.1 rated frames
- Innovative design for Unparalleled Protection
- Designed for Superior Protection & Performance
- Weight: 76 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Matte Black



BRICK



- Peak Performance Ballistic Eyewear
- ANSI Z87.1 rated frames
- Innovative design for Unparalleled Protection
- Designed for Superior Protection & Performance
- Weight: 56 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Matte Black, Crystal Metallic



ENZO



- Peak Performance Ballistic Eyewear
- ANSI Z87.1 rated frames
- Innovative design for Unparalleled Protection
- Designed for Superior Protection & Performance
- Weight: 66 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Matte Black

NIKE Eyewear Collection

BOOST



- Fabricated from lightweight polycarbonate
- Soft rubber nose pads and temple tips for stability and grip
- Fits medium to large faces
- Weight: 74 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Anthracite



FUEL

- Bio-injected frame material contains at least 40% castor bean oil
- Soft rubber nose pads and temple tips
- Fits small to medium faces
- Weight: 56 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Black



PREMIER

- Lightweight nylon frame for comfort, durability
- 6-base frame design for maximum coverage
- Fits medium to large faces
- Weight: 78 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Black, Tortoise



RABID

- Lightweight frame design provides both maximum coverage & clarity
- Secure wrap temple arms for stability, comfort
- Fits medium to large faces
- Weight: 64 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Black

All Rudy Project frames are entirely designed, engineered, and made in Italy using advanced materials. Rudy Project frames blend cutting-edge ergonomic features with advanced aero spatial materials such as proprietary Kynetium, Grilamid, Titanium & Carbon Fiber to create extremely light-weight versatile and high performance glasses.



CONTENDER



- Lightweight Titanium & Carbon frame
- Ergonose System for superior fit
- Integrated Side Shields for maximum protection
- Weight: 73 Grams
- Front .75mm Pb Side .5mm Pb Protection
- Frame Colors: Black



VOLTAGE



- Ultralight Titanium Frame
- Ergonose System for maximum comfort
- Weight: 59 Grams
- Front .75mm Pb Protection
- Frame Colors: Carbon Grey



SENSOR



- Ultralight Hybrid Titanium & Kynetium frame
- Ergonose System for extraordinary fit
- Ideal for all Rx types
- Weight: 54 Grams
- Front .75mm Pb Protection
- Frame Colors: Carbon Silver

NEW PRODUCT

PRIVATE LABEL COLLECTION

BARRIER SIREN



GLOSS BLACK

Side Shield Protection



- Lightweight, durable frame design
- Ideal for small faces
- Ventilated nose bridge reduces fogging
- Integrated side shields for enhanced lateral protection
- Weight: 65 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Colors: Gloss Black, Polished Clear



POLISHED CLEAR



Available Rx's

PL SV BF PG

NEW PRODUCT

PRIVATE LABEL COLLECTION

BARRIER BRAZEN



PL-BRAZEN-T

TORTOISE



BLACK



PL-BRAZEN-B



Side Shield Protection



- Universally fitting frame for all face sizes
- Lightweight and stylish design
- Available with side shields
- Ventilated nose bridge for additional comfort
- Weight: 76 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Colors: Black, Tortoise

Available Rx's





VELOCITY



- The unique materials & geometries of this goggle afford the wearer superior splash protection, heat & impact resistance.
- Clear Zone Lens treatments provide an optimal combination of anti-scratch coating on the outside & anti-fog coating on the inside
- Low profile frame maximizes peripheral vision
- Full-perimeter ventilation and filtration minimizes lens fogging



DET CORD



- Unobtanium® nose piece & ear socks for enhanced grip
- Three-Point Fit holds lenses in precise optical alignment while eliminating pressure points
- Weight: 74 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Colors: Matte Black



GASCAN



- O-Matter™ stress-resistant frame material is both lightweight & durable for all-day comfort and protection
- Three-Point Fit holds lenses in precise optical alignment while eliminating pressure points
- Weight: 74 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Colors: Polished Black, Matte Black, Multi-Camo



ASTRO DLX

- Enhanced strength and durability
- Wire-core reinforced temple arms
- Flexible nose pads design to fit any facial profile
- Available with side shields
- Weight: 67 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Color: Blue/Black, Black/Black, Red/Black



BLADE II

- Adjustable nose piece
- Integrated side shields for enhanced lateral protection
- Rx Capable
- Weight: 52 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Black



SLIM-T

- Precision machined ultralight frame
- Dual hinges for ultimate comfort and flex
- Designed for all face types
- Integrated .5mm Pb side shields
- Weight: 66 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Gun Metal Grey





MX-30

- Features a super soft fixed rubber nose piece
- Weight: 63 Grams
- Lead equivalent Front Lens: .75 mmPb
- Optional goggle style w/ adj. head strap
- Available with side shields
- Frame Colors: Crimson, Black, Silver, Blue



NOW AVAILABLE IN GLOBAL FIT
Frame Color: Black with Red Accents



FITOVER

- Perfect to fit over existing personal Rx eyewear
- Goggle design protects from fluid splashes
- Weight: 83 Grams
- Lead equivalent Front Lens: .75 mmPb
- Lead equivalent Side Lens: .5 mmPb:
- Frame Colors: Black, Tortoise



FITOVER XL

- Designed for larger faces and glasses
- Goggle design protects from fluid splashes
- Weight: 95 Grams
- Lead equivalent Front Lens: .75 mmPb
- Lead equivalent Side Lens: .5 mmPb:
- Frame Colors: Black, Tortoise



BARRIER LITE™

- Polycarbonate wrap around frame
- Weight: 60 Grams
- Lead equivalent Front Lens: .75 mmPb
- Available with side shields
- Frame Colors: Crimson, Black, Silver



T-ZONE

- Available in 3 sizes; small, medium, large
- Fits all facial types from small to X-large
- Fashionable wrap around style
- Optional goggle style with adj. head strap
- Weight: 66 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Colors: Black, Blue, Red, Silver



53 WRAP

- Flexible hinges
- Weight: 87 Grams
- Lead equivalent Front Lens: .75 mmPb
- Lead equivalent Side Lens: .5 mmPb
- Frame Colors: Black Crystal, Tortoise



BLADE II

- Adjustable nose piece
- Integrated side shields for enhanced lateral protection
- Rx Capable
- Weight: 52 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Black



CONTENDER



- Lightweight Titanium & Carbon frame
- Ergonose System for superior fit
- Integrated Side Shields for maximum protection
- Weight: 73 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Black, Graphite Red



SENSOR



- Ultralight Hybrid Titanium & Kynetium frame
- Ergonose System for extraordinary fit
- Ideal for all Rx types
- Weight: 54 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Carbon Silver



VOLTAGE



- Ultralight Titanium Frame
- Ergonose System for maximum comfort
- Weight: 59 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Carbon Grey



PL SV BF PG

CROSSLINK XL

- Sleek and modern look matches Oakley performance innovation with a comfortable fit
- Constructed of light weight, durable O-Matter[™] material
- Unobtainium[™] nosepads and ear socks for enhanced grip
- Weight: 62 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Colors: Grey, Satin Black



PL

MAINLINK

- Lightweight O-Matter[®] frame for optimal comfort
- Three-Point fit holds the lenses for precise optical alignment
- Integrated Unobtainium for secure grip and fit
- Weight: 66 Grams
- Lead Equivalent Front Lens: .75mm Pb
- Frame Colors: Matte Black



PL SV BF PG

TWOFACE

- Unique pairing of lightweight O-Matter[™] upper & high strength Aluminum lower rims
- Lens positioning offers superior comfort, clarity, & razor-sharp vision at every angle
- O-Matter[™] stress-resistant frame material is lightweight & durable for comfort and protection
- Weight: 64 Grams
- Lead equivalent Front Lens: .75 mmPb
- Frame Colors: Steel w/Grey, Polished Black



PL SV BF PG

MX-30

- Features a super soft fixed rubber nose piece
- Weight: 63 Grams
- Lead equivalent Front Lens: .75 mmPb
- Optional goggle style w/ adj. head strap
- Available with side shields
- Frame Colors: Black & Red

ClearView AR™

The Ultimate in Clarity



Fog Free Anti-reflective Coating




How does **ClearView AR™** work?

As light enters a lens, roughly 85% passes through while the remaining light is reflected off the surface as glare. Glare can be distracting, especially in a healthcare environment where lighting is intense. Realview AR coating will allow 99.5% of light to pass through resulting in superior clarity.

ClearView AR™ Benefits

- Improves Your Vision
- Scratch Resistant
- Reduces Fingerprints and Smudging
- Repels Oils
- Easier to Clean



Leaded Eyewear

PLEASE NOTE:

In excess of the following powers cannot be manufactured

+/- 8.00 Rx maximum power.





























+/- 5.00 Rx power must use a lenticular design.

+/- 4.00 Rx power must include ClearView AR[™] anti-reflective coating.

+/- 4.00 Rx maximum power for all Bifocals.

Rx Limitations

 SV	Single Vision
 BF	Bifocal
 PG	Progressive

Frame	Rx Limitations SV or PG	Limitations BF	Available Lens
Astro DLX	+/- 4.00 Total Power	+/- 4.00 Total Power	
Atomic	+/- 6.00 Total Power	+/- 4.00 Total Power	
Barrier Lite	+/- 2.00 Total Power	N/A	
Blade II	+/- 6.00 Total Power	+/- 4.00 Total Power	
Brazen 	+/- 3.00 Total Power	+/- 3.00 Total Power	
Brine 	+/- 2.00 Total Power	+/- 2.00 Total Power	
Circuit XL	+/- 8.00 Total Power	+/- 4.00 Total Power	
Crosslink 	+/- 3.00 Total Power	+/- 3.00 Total Power	
Grid II (minimum PD: 31)	+/- 8.00 Total Power	+/- 4.00 Total Power	
Holbrook Rx 	+/- 8.00 Total Power	+/- 4.00 Total Power	
Jupiter Squared 	+/- 4.00 Total Power	+/- 3.00 Total Power	
Justin 	+/- 3.00 Total Power	+/- 3.00 Total Power	
MX-30	+/- 8.00 Total Power	+/- 4.00 Total Power	
New Wayfarer 	+/- 4.00 Total Power	+/- 3.00 Total Power	
Sensor 	+/- 8.00 Total Power	+/- 4.00 Total Power	
TwoFace 	+/- 3.00 Total Power	+/- 3.00 Total Power	
T-Zone	+/- 8.00 Total Power	+/- 4.00 Total Power	
Voltage 	+/- 2.00 Total Power	N/A	

See All of Our Radiation Protection Products at:
WWW.BARRIERTECHNOLOGIES.COM



Aprons



Leaded Eyewear



Gloves



Thyroid Collars



Mobile Barriers



Accessories



Secure Shields



Dental



7060 W. STATE ROAD 84, Suite 8, DAVIE, FL 33317

Tel: (800) 353-4350 ▪ Fax: (800) 353-4351

Tel: (954) 921-0330 ▪ Fax: (954) 921-0331

WWW.BARRIERTECHNOLOGIES.COM