

PAXIS 250 AUTO



*The Paxis 250 is a universal radiographic solution
All positions can be called up automatically
or done manually through motorization functions*

PAXIS 250 AUTO

The Paxis 250 is one of the most versatile X-Ray positioning available

Having all the conveniences of a universal positioning device combined with ease of use of three preprogrammed SID positions (40" and 72" in the horizontal and 40" in the Vertical), the Paxis 250 is an ideal choice for any imaging room.

Standard configuration is full motorized positioning and auto-positioning is built in

This means that you now have a system with tremendous ease of operability

Our newly added auto-positioning feature streamlines operator performance

It has a motorized source to image distances of 40" to 72"

Additional motorization is found in the vertical travel as well as the arm rotation



***The Paxis 250 Auto-Positioning
brings unrivaled performance in
a compact package that is ideal
for all users and needs***

***Its full range of motions
means you can perform
skull procedures all the way
through to the
lower extremities***

***The Paxis 250 Auto has
a tremendous amount
of flexibility to suit your
every need***



One solution for all your needs

Accogent Pausch • 808 Shrewsbury Avenue • Tinton Fall, NJ • T: 732.747.6110 • accogent.com