Radmedix

Radmedix 2510 Lance Drive Dayton, OH 45409
Tel: 844.723.6334 | Fax: 513.278.0150 | info@radmedix.com
© 2025 Radmedix Inc. All rights reserved.



Radmedix

ACUITY DR GO

Equine & Mixed Practice Solution

The Gold Standard of Veterinary X-Ray Imaging.

Acuity DR Go — Portable Veterinary DR System

The *Acuity DR Go* veterinary solution has everything you want from a versatile, portable DR system. Specifically designed for equine and mixed practice imaging, it's easy to transport in the self powered, portable case without worrying about power cords to operate. It can be configured to any clinical needs with either the Acuity DR 10x13 or 14x17 detector options.

KEY FEATURES:

- Lightweight Durable Hard Shell Case with Carrying Strap
- Built-In PC Workstation with Touchscreen Monitor
- Self-Powered by Dual Dewalt Batteries
- Wireless Control from Smartphone App
- Intuitive AccuVue Vet Software
- High Definition Imaging



Lightweight High-Tech Composite

Extremely durable lightweight case made with the latest high-tech composite material with shock absorbing edge reinforcement, also splash proof (IPX4).



Built-In X-Ray PC Workstation

Built-in high performance PC workstation with Windows operating system. Also comes with a built-in keyboard and trackpad mouse.

Designed for Equine and Mixed Practices.

Huge 22” HD Touchscreen Monitor

The Acuity DR Go features a large, anti-glare tilting 22" touchscreen display when taking images on the go as well as a built-in trackpad mouse for precise post processing tools.



Padded DR Compartment

Easy removal and storage of the x-ray detector. Fits panel and protection case.

Power Saving Mode

Fold down the monitor to activate power saving mode and continue where you left off within a set timeframe. System auto shuts down if not reopened.

USB Slot





Intelligent Battery Design

The Acuity DR Go has a unique battery docking system that allows you to swap in power tool rechargeable batteries on the go, to keep the case ready for clinical use for hours without recharging. *Page 3 for more details

All Day Clinical Use with Intelligent Battery Design.

Intelligent Battery Design - Battery Swap

The Acuity DR Go case uses a completely new concept for the power supply, it allows you to use standard rechargeable batteries to keep the case powered indefinitely.



ADVANTAGES:

- ✔ *Batteries can be changed during operation. No shutdown of the system is necessary.*
- ✔ *Batteries are not discharged when the case is switched off. Case is ready for immediate use, even if it has not been used for a long time.*
- ✔ *Unit automatically switches to stand-by mode when the lid is closed. After changing location, the computer does not have to be restarted.*
- ✔ *The standard rechargeable batteries can be bought in almost any DIY store worldwide.*





Acuity DR Go Comes with 2 Dewalt Rechargeable Batteries

(2) 20V/5 Ah for a operating time of approx. 8.5 hrs for 500 x-ray exposures

Configure Acuity DR Go to Your Exact Clinical Needs.

PAIR WITH OUR ACUITY DR 10x13



Our *Acuity DR 10x13 panel* is our "work horse" in superior equine imaging. Perfect for smaller format detailed exams.

Features: *100 micron pixel pitch for sharper image detail, extended battery life of up to 8 hrs, and IPX5 rated for resistance against liquid intrusion.*

PAIR WITH OUR ACUITY DR 14x17



Our *Acuity DR 14x17 panel* is sold into the most demanding veterinary environments requiring unique imaging versatility.

Features: *100 micron pixel pitch for sharper image detail, extended battery life of up to 9.5 hrs, and IP56 rated for resistance against liquid and dust intrusion.*

Carrying Strap Included for Comfortable Transportation



Intuitive Operation with AccuVue Vet Software.

Acuity DR Go — AccuVue Vet Software

Radmedix's powerful AccuVue Vet software offers a clean, customizable user interface for fluid workflow even when complex exams are demanded. It's 10 point touch optimization offers the look and feel of a modern software GUI users expect from a seasoned imaging system. AccuVue Vet was designed with premium equine and mixed practice veterinary imaging in mind.



KEY FEATURES:

- Modern graphical user interface (GUI) available in almost any language with touchscreen operation to ensure a quick, efficient, and smooth workflow
- Use of DICOM Procedure Codes for the transfer of all relevant examination data directly from the connected patient management system
- Safe and fast registration of emergency patients allowing user to switch between projections to avoid having to re-position the patient frequently
- Capture of patient data via DICOM Worklist, BDT/GDT, HL7 or other protocols – data may also be captured manually
- Freely configurable body parts with more than 200 projections and numerous possible adjustments already included
- Allows the user to subsequently add images to an examination, even after that examination has been completed + more

Acuity DR Go Specifications | What's Included?



RadmediX

Acuity DR Go — Specifications

▪ Extremely durable case with edge shock absorption, splash-proof (IPX4), and is lockable	▪ Comes with AccuVue Vet software
▪ Anti-glare 22" full HD touchscreen monitor, keyboard, and trackpad mouse	▪ Configure with either our 10x13 or 14x17 detector, includes DR protection case
▪ Carrying strap for comfortable transport	▪ (2) 20V/5 Ah Dewalt rechargeable batteries
▪ Padded DR compartment for 10x13 or 14x17 detector	▪ Use with your smartphone for wireless control
▪ Connections: 1 x USB, 1 x LAN, 1 x WLAN, 1 x Bluetooth, 1 x DC socket	▪ High performance PC with Windows OS
▪ Total Weight: 21 lbs (includes complete x-ray case, monitor, and detector)	▪ Dimensions: 21" W x 24" H x 20" D